

## What Are Medical Gases?

### Oxygen

Oxygen may be used for patients requiring supplemental oxygen via a mask. Usually accomplished by a large storage system of liquid oxygen at the hospital which is evaporated into a concentrated oxygen supply, pressures are usually around 380 kPa (55 psi). In small medical centers with a low patient capacity, oxygen is usually supplied by multiple standard cylinders.

Medical oxygen is used:

For virtually all modern anesthetic techniques as well as pre-and post-operative management. Oxygen provides life support by restoring tissue oxygen levels in a range of conditions, including:

- Cyanosis as a result of cardio-pulmonary disease
- Surgical trauma, chest wounds and rib fractures
- Shock, severe hemorrhage and coronary occlusion
- Carbon monoxide poisoning
- Hyperpyrexia
- Major trauma, such as road accidents and gunshot wounds.

Oxygen plays a vital role in the management of sudden cardiac and respiratory arrest - whether drug induced or traumatic - and in the resuscitation of critically ill patients when circulation is impaired. It is also used in neo-natal resuscitation.

### Nitrous Oxide

Nitrous Oxide is supplied to various surgical suites for its anesthetic functions during pre-operative procedures. Delivered to the hospital in standard tanks and supplied through the Medical Gas system. System pressures around 345 kPa (50 psi).

Medical nitrous oxide is used:

- For the relief of severe pain, usually in emergency situations, by inhalation with 50% oxygen
- During induction and maintenance of anesthesia, in controlled situations

### Nitrous Oxide and Oxygen

A mix of nitrous oxide 50% and oxygen 50% is a medical anesthesia gas, commonly known as Entonox (a registered trademark of BOC) or Nitronox, or colloquially as “gas and air”, and is frequently used in pre-hospital care, childbirth and emergency medicine situations by medical professionals such as doctors, nurses, midwives and paramedics.



VALVES

CYLINDERS

CRYOGENIC

HOSES &amp; PAINTS

ACCESSORIES

Evergreen Midwest Co  
8976 Osborne Drive  
Mentor, OH 44060  
Toll: 800.659.3358  
Phone: 440.255.5540  
Fax: 440.255.6434  
Email: sales@emwco.com

[www.EverGreenMidwest.com](http://www.EverGreenMidwest.com)

### *Nitrogen*

Nitrogen is typically used to power surgical equipment during various procedures. Pressures range around 1.2 MPa (175 psi) to the various locations.

### *Carbon Dioxide*

Typically used to inflate or suspend tissues during surgery, and also used in laser surgeries. System pressures are maintained at about 345 kPa (50 psi).

Medical carbon dioxide is used:

- To increase the depth of anesthesia rapidly when volatile agents are being administered. It increases depth of respiration and helps to overcome breath-holding and bronchial spasm
- To facilitate blind intubation in anesthetic practice
- To facilitate vasodilation, and thus lessen the degree of metabolic acidosis during induction of hypothermia
- To increase the cerebral blood flow in arteriosclerotic patients undergoing surgery
- To stimulate respiration after a period of apnea
- In chronic respiratory obstruction after it has been relieved
- To prevent hypercapnia during hyperventilation
- For clinical and physiological investigation.

### *Medical Air*

Medical Air is supplied by a special air compressor to patient care areas using clean outside air. Pressures are maintained around 380 kPa (55 psi). Medical air can be used as surgical air when the pressures are raised to 7 Bars. Surgical air is majorly used in driving pneumatic tools like dental hand pieces. Medical air should never be used as oxygen.

Medical air is used:

- As a replacement for atmospheric air when the atmosphere is contaminated by noxious fumes, vapors or gases
- In anesthesia as a carrier for volatile anesthetic agents
- As a power-source for pneumatic equipment
- In ventilators and incubators to provide uncontaminated and controlled air flow.

### *Medical Test Gas Mixtures*

There are many gas mixtures used for clinical and medical applications. They are often used for patient diagnostics such as lung function testing or blood gas analysis. Test gases are also used to calibrate and maintain medical devices used for the delivery of anesthetic gases.

### *Medical Culture Growth Mixtures*

Culture growth applications include controlled aerobic or anaerobic incubator atmospheres for biological cell culture or tissue growth. Controlled aerobic conditions are created using mixtures rich in oxygen and anaerobic conditions are created using mixtures rich in hydrogen or carbon dioxide.

### *Quick Facts:*

#### **Air (AIR)**

Colorless, odorless, tasteless. Chemical composition is 78% nitrogen, 21% oxygen, and 1% trace elements.

#### **Oxygen (O<sub>2</sub>)**

Colorless, odorless, and tasteless medical gas. Not flammable but does support combustion

#### **Carbon Dioxide (CO<sub>2</sub>)**

Colorless, odorless and tasteless medical gas. Does not support combustion or life.

#### **Helium (He)**

Colorless, odorless and tasteless medical gas. Normally used as heliox in a 80/20 or a 70/30 mixture

#### **Nitric Oxide (NO)**

Colorless and has a faint taste and smell. Used as a vasodilator and used more in the neonatal unit

#### **Nitrous Oxide (N<sub>2</sub>O)**

Colorless, odorless and tasteless medical gas. Used as an anesthetic

#### **Ways to Manufacture O<sub>2</sub>**

Fractional Distillation, Electrolysis of H<sub>2</sub>O, Chemical Decomposition, Physical Separation

#### **Fractional Distillation**

Cheapest and most common method to manufacture O<sub>2</sub>. Air is filtered to remove pollutants, H<sub>2</sub>O and CO<sub>2</sub>. Air is then liquefied by compression and cooled by rapid expansion. The resulting mixture is heated in a distillation tower and then transferred to cryogenic storage cylinders

#### **Electrolysis of H<sub>2</sub>O**

Separates oxygen from water



VALVES

CYLINDERS

CRYOGENIC

HOSES & PAINTS

ACCESSORIES

Evergreen Midwest Co  
8976 Osborne Drive  
Mentor, OH 44060  
Toll: 800.659.3358  
Phone: 440.255.5540  
Fax: 440.255.6434  
Email: [sales@emwco.com](mailto:sales@emwco.com)

[www.EverGreenMidwest.com](http://www.EverGreenMidwest.com)

### Chemical Decomposition

Heating of oxygen to remove impurities

### Physical Separation

Uses oxygen concentrators to extract oxygen from ambient air; produces the least amount of oxygen

### Most Frequently Used Cylinders in Care Settings

H, G, and E cylinders

### Color Coding for Cylinders

<b>Air</b> Yellow or black/white	<b>Oxygen</b> Green or white *green used most frequently*	<b>Carbon Dioxide</b> Gray	<b>Carbon Dioxide/Oxygen</b> Gray/green
<b>Helium</b> Brown	<b>Nitrogen</b> Black	<b>Nitrous Oxide</b> Light blue	

### Gas Safety Systems

PISS

DISS

\*not universal systems\*

**Pin Index Safety System** - PISS

**Diameter Index Safety System** - DISS

### Types of Pressure Relief Devices Used for Gas Cylinders

Rupture disk (frangible disk)

Fusible plug

Spring loaded

### Medical Gas Distribution Systems

Bulk systems

Portable cylinders

Liquid systems

**Central Supply Systems**

Continuous  
Alternating

**Continuous Supply System**

2 or 3 tanks that are consistently running  
Alternating Supply System

1 backup tank in place at the care setting used in case the continuous supply systems fail

**Working Gas Pressure**

50 pounds per square inch gauge (psig)

**Piping Systems**

Gas delivery system from the central supply units to sites or zones within the hospital; allows for the quick isolation of all independent areas in case of fire with the shut off valves

**Determining Duration for Flow for a Gas Cylinder**

$CF \times PSIG \div LPM$

cf= cylinder factor

psig= pounds per square inch gauge

lpm= liters per minute

**Medical Air Supply**

Piston air compressors  
Diaphragms  
Rotaries

**Reducing Valve**

Device that reduces high pressure gases from cylinders or bulk storage units to lower the working pressure

**Flowmeters**

Device used to control and indicate gas flow, usually in LPM

Thorpe Type  
Bourdon  
Flow Restrictor

**VALVES****CYLINDERS****CRYOGENIC****HOSES & PAINTS****ACCESSORIES**

Evergreen Midwest Co  
8976 Osborne Drive  
Mentor, OH 44060  
Toll: 800.659.3358  
Phone: 440.255.5540  
Fax: 440.255.6434  
Email: [sales@emwco.com](mailto:sales@emwco.com)

[www.EverGreenMidwest.com](http://www.EverGreenMidwest.com)

**Thorpe Tube**

Most commonly used. Uses a float and a needle valve to control the flow of a gas

**Bourdon Flowmeter**

A reducing valve that controls the pressure gradient across an outlet with a fixed orifice. So long as the pressure distal to or downstream from the orifice remains atmospheric, the indicated flow is accurate

**Flow Restrictor**

Operates on the same principles as the Bourdon flowmeter. Gas flow can be increased by raising driving pressure across a fixed resistance

\*used for lpm of  $\frac{1}{2}$ -3\*

**Oxygen Adder**

Consists of 2 flowmeters; the flow is directed through a humidifier and then to the patient. The  $\text{FiO}_2$  rate is dependent on the air entrainment ratio

\*simplest method\*

**Oxygen Blender**

Uses working pressure source (50 psig) of air and oxygen to deliver precise  $\text{FiO}_2$ . Works pneumatically

\*compact and convenient\*

**Components of an Oxygen Blender**

Alarm module

Pressure-balancing module

Proportioning module

**Problems with Oxygen Blenders**

Contaminated gas sources,  $\text{H}_2\text{O}$ , particulate matter, corrosion (prevents check valves from seating properly and diaphragms from properly moving), uses filters at inlets

DISCLAIMER: The information contained herein is provided as a service to the public and is intended for informational purposes only. Evergreen Sales and Services Corp. expressly warns that information contained herein may not be completely accurate and/or up to date. Therefore, neither Evergreen Sales and Services Corp. nor any of its affiliates assumes any liability whatsoever for the accuracy or completeness of the information contained herein. The statements herein shall have no force or effect. If you detect any omissions, misstatements or errors, please contact us immediately.



VALVES

CYLINDERS

CRYOGENIC

HOSES &amp; PAINTS

ACCESSORIES

Evergreen Midwest Co  
8976 Osborne Drive  
Mentor, OH 44060  
Toll: 800.659.3358  
Phone: 440.255.5540  
Fax: 440.255.6434  
Email: sales@emwco.com

[www.EverGreenMidwest.com](http://www.EverGreenMidwest.com)